COVIDIEN™

MAHURKAR™ Elite

Acute Triple Lumen Catheter

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Sterilized with ethylene oxide. Sterile and non-pyrogenic in unopened and undamaged package.

DO NOT RESTERILIZE.

CAUTION: This device is for single use only. Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

READ ALL INSTRUCTIONS, WARNINGS, AND CAUTIONS CAREFULLY PRIOR TO USE.

DESCRIPTION

The MAHURKAR Elite catheter is a radiopaque, polyurethane tube with two clear polyurethane catheter extensions and one clear polyurethane infusion lumen. The three internal lumina can be distinguished by the color-coded luer-lock adapters on the extensions:

red adapter = arterial lumen

blue adapter = venous lumen

white adapter = medial lumen

The arterial lumen provides outflow from the patient; the venous lumen provides return; the medial lumen is for infusion of fluid, blood products, medications, blood sampling, pressure injection of contrast media and central venous pressure monitoring.

The catheter comes in a variety of sizes which are listed below. Options include curved or straight extensions. Refer to component specific instructions in the following pages, component packaging, and/or the product insert (if provided).

INFUSION FLOW RATES		
Catheter Length	Infusion Pressure (cm of water)	Medial Lumen, White (mL per min)
13 cm	100	72.1
16 cm	100	63.6
20 cm	100	59.7
24 cm	100	56.6
30 cm	100	49.0

DYNAMIC FLOW RATES

Flow Rate (mL/min)	Lumen	Average Pressure (mmHg)				
		13 cm	16 cm	20 cm	24 cm	30 cm
100	Venous	6	6	6	8	9
	Arterial	-5	-6	-7	-8	-10
150	Venous	10	10	12	14	16
	Arterial	-10	-12	-13	-14	-17
200	Venous	15	16	19	21	25
	Arterial	-17	-19	-22	-25	-30
250	Venous	21	24	29	32	39
	Arterial	-25	-29	-33	-37	-46
300	Venous	29	34	41	47	56
	Arterial	-34	-40	-47	-52	-64
350	Venous	38	44	53	62	74
	Arterial	-45	-53	-61	-68	-84
400	Venous	49	55	68	78	93
	Arterial	-58	-67	-77	-87	-108
450	Venous	61	68	82	96	113
	Arterial	-72	-83	-95	-107	-133
500	Venous	74	81	99	115	135
	Arterial	-86	-99	-115	-129	-161

INDICATIONS

The MAHURKAR Elite triple lumen catheter is intended for short-term central venous access for hemodialysis, apheresis, infusion, central venous pressure monitoring and pressure injection of contrast media. The maximum recommended infusion rate is 5 mL/sec for power injection of contrast media.

CONTRAINDICATIONS

The catheter is not intended for any purpose other than indicated in these instructions. The catheter should not be implanted femorally longer than three to four days maximum per K-DOQI guidelines. Do not use this catheter in thrombosed vessels or for subclavian puncture when a ventilator is in use.

POTENTIAL COMPLICATIONS

Subclavian, jugular, femoral, or saphenous vein insertion: hematoma; central venous thrombosis/stenosis; sepsis; catheter thrombosis; exit site infection; trauma to major vessel or right atrium; air embolism; cardiac tamponade; hemorrhage; hemothorax; pneumothorax; cardiac arrhythmia; carotid sinus tract infection; brachial plexus injury.

Femoral or saphenous vein insertion: retroperitoneal bleed

CLINICAL SUMMARY

A randomized, controlled, multi-centered clinical study was conducted to compare the safety and efficacy of Mahurkar Triple Lumen Catheter, 12 Fr. to the Mahurkar Dual Lumen Catheter, 11.5 Fr in patients undergoing acute hemodialysis or apheresis procedures. The primary end point of this study was catheter related blood stream infections. Device safety was determined by the incidence of complications, specifically catheter-related infections, occurring with the use of the MAHURKAR Triple Lumen Catheter. Device effectiveness was determined by blood flow performance during hemodialysis or apheresis treatment. Patients were observed for a maximum of 30 days, which included a maximum catheter insertion period of 29 days and a minimum observation period of 24 hours following catheter removal.

In this study, subjects were included for the following criteria: met labeled indication for hemodialysis or apheresis, had a life expectancy of at least 30 days, were able to give written informed consent and weighed ≥ 40 kg. The subjects were excluded based on the following criteria: required ventilators under positive end expiratory pressure mode at values > 5 and required catheter placement in the subclavian or jugular vein, required an exchange from an existing non study catheter for the study catheter to support eligibility, were previously enrolled in the MAHURKAR Triple Lumen Catheter clinical study. Four hundred and eighty-five patients (485) were enrolled in the study from 11 sites. The following table lists the patient demographics information.

		All Patients (n = 485)	Dual Catheter (n = 243)	Triple Catheter (n = 242)
Gender				
	Female	240 49%	116 48%	124 51%
	Male	245 51%	127 52%	118 49%
Age*				
	Median (range)	59 (18-95)	61.5 (21-92)	56.5 (18-95)
		84 17%	34 14%	50 21%
		161 33%	79 33%	82 34%
		237 49%	127 53%	110 45%
Hx Diabetes				
	No	271 56%	135 56%	136 56%
	Yes	214 44%	108 44%	106 44%
Immunosuppressant Use				
	No	443 91%	218 90%	225 93%
	Yes	42 9%	25 10%	17 7%
Other Risk Factors (surgical incisions, skin trauma, other)				
	No	337 69%	170 70%	167 69%
	Yes	148 31%	73 30%	75 31%
Number of Invasive Devices Pre-Catheter Insertion				
	0	141 29%	64 26%	77 32%
	1	165 34%	84 35%	81 33%
	2 or more	179 37%	95 39%	84 35%
Pre-Insertion Infection				
	Yes	89 18%	46 19%	43 18%
	No	271 56%	134 55%	137 57%
	Unknown	125 26%	63 26%	62 26%
Patient Incontinence**				
	No	205 43%	96 41%	109 46%
	Yes	268 57%	139 59%	129 54%
*Age was missing for 3 dual lumen catheter patients				
**N=473 applicable to femoral insertions only.				

In summary, out of the 393 evaluable patients studied, the catheter related blood stream infection rate was 7.5% in the MAHURKAR Triple Lumen Catheter group, in comparison to the rate of 8.3% in the MAHURKAR Dual Lumen Catheter group. The difference was not statistically significant (P=0.77 by Chi-square Test).

WARNINGS AND PRECAUTIONS:

- The catheter should be inserted and removed only by a qualified, licensed physician or other health care practitioner authorized by and under the direction of such physician.
- The medical techniques and procedures described in these instructions do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.
- Any complications arising from non-hemodialysis/apheresis procedural use of the catheters are the sole responsibility of the person(s) using the catheter.
- Do not use the catheter if package has been damaged or previously opened. Do not use the catheter if it is crushed, cracked, cut, or otherwise damaged.
- The catheter should not remain in place for more than 29 days.
- Discard after use: the catheter is for one time insertion only.
- This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.
- Observe and follow all expiration dates for all individually packaged components.
- Observe sterile technique at all times when handling the catheter.
- Use the guidewire straightener to insert the "J" end of the guidewire into the introducer needle. **Do not insert or withdraw the guidewire forcibly from any component; the wire could break or unravel.**
- **Immediately after insertion and before using the catheter, use X-ray or fluoroscopy to verify that the catheter tip is positioned correctly. For jugular and subclavian insertions, ensure that the tip is positioned at the junction of the superior vena cava and the right atrium. Failure to do so could result in serious trauma or fatal complications.**
- To avoid exposure to blood borne pathogens, use universal precautions with infection control and sharp safety components.
- Do not infuse incompatible drugs simultaneously through the same lumen; precipitation could occur.
- Use of the subclavian vein for catheter placement may result in subclavian vein stenosis. Subclavian vein stenosis may impair use of the ipsilateral extremity for future vascular access. Use of the jugular vein may be preferable.
- Avoid air embolism by keeping the catheter extension tubing clamped at all times when not in use and by filling the catheter with sterile saline prior to implantation. With each tubing change purge air from the tubing and aspirate air from the catheter.
- Overtightening the catheter connections can crack the adapters.
- Do not clamp the triple lumen tubing portion of the catheter; clamp only the extensions. When clamping, use only clamps supplied with the catheter.
- Clamping the catheter extensions repeatedly in the same spot could weaken the tubing: change the position of the clamp regularly to prolong the life of the tubing. Avoid clamping near the adapter and hub.
- Exercise caution when using sharp instruments near the catheter. The catheter tubing can tear when subjected to nicks, excessive force, or rough edges.
- Inspect the catheter frequently for nicks, scrapes, cuts, etc. which could impair performance.
- Do not use alcohol or any solution containing alcohol to lock the catheter between dialysis treatments.
- When instilling heparin to lock the catheter, flush quickly and clamp immediately to ensure heparin reaches the distal end of the lumen. Do not infuse against a closed clamp or forcibly infuse a blocked catheter.

- Per KDOQI guidelines, remove the catheter when no longer needed. It is recommended that internal jugular catheters should be used no more than one week and femoral catheters should be used no more than five days in bed-bound patients with good exit site care.
- Do not use acetone or any solution containing acetone on any part of the catheter. Exposure to these agents may cause catheter damage.
- Any ONE of the following can be used as a cleaning agent: aqueous-based povidone-iodine, chlorhexidine, amukin 50%, bacitracin ointment, hydrogen peroxide, Betadine™, Exsept™, 70% isopropyl alcohol, Chloraprep™, Neosporin™, Octenisept™, Dilutis™, Polysporin™, Bactroban™ cream, aqueous chlorhexidine gluconate or chlorhexidine 2% with alcohol. PLEASE NOTE: mixing of these solutions has not been tested and is not recommended.
- It is recommended that only luer-lock (threaded) connections be used (including syringes, bloodlines, IV tubing, and sealing caps) with the catheter.
- To prevent damage to the catheter and ensure patency before high pressure injection, forcefully flush the catheter using normal saline with a 10 mL or larger syringe to clear any potential occlusion. Resistance to flushing may indicate partial or complete catheter occlusion. Do not proceed with high pressure injection procedure until occlusion has been cleared.
- Prior to high pressure injection, failure to ensure patency of the catheter may result in catheter malfunction.
- Use only the medial lumen with the white luer adapter for high pressure injection of contrast media. Do not use arterial (red) or venous (blue) lumens for high pressure injection.
- With an occluded catheter the pressure limiting feature for a power injector machine may not prevent over-pressurization of the catheter, which may lead to catheter failure.
- There is potential for catheter failure or catheter tip displacement if the maximum flow rate of 5 mL/sec. is exceeded.
- Not all patients are appropriate for high pressure injection even though the catheter can withstand the procedure. It is necessary to have a trained clinician responsible for assessing the health status of a patient that is considered a candidate for a high pressure injection procedure.
- Potential catheter failure may result from high pressure injection of a catheter that has signs of clavicle-first rib compression or pinch-off.
- High pressure injection should be stopped immediately if signs of swelling, local pain, or signs of extravasation are noted.

INSERTING THE CATHETER

Supplies Needed

- | | |
|--|--|
| (1) catheter | (1) fenestrated sterile drape |
| (1) 5 mL ampule lidocaine 1% | (1) mask |
| (1) #11 scalpel | (1) gown |
| (1) suture with curved needle | (1) pair gloves |
| (1) needle holder | (1) cap |
| (3) needles -18 G (1.2 mm); 22 G (0.7 mm); 25 G (0.5 mm) | (2) wound dressings |
| (2) 5, 6, or 10 mL syringes | (1) razor (optional) |
| (3) 3 mL syringes | (1) pkg. appropriate antiseptic |
| (1) 18 G (1.2 mm) introducer needle | (4) 4" x 4" (10.2 cm x 10.2 cm) gauze sponges |
| (3) sealing caps | (1) 0.035" (0.889 mm) J/straight stainless steel guidewire |
| Appropriate sized dilators | |

Heparinized saline: Use concentrations approved by your institution

PROCEDURE

CAUTION: READ THIS PROCEDURE CAREFULLY BEFORE INSERTING THE CATHETER.

The operating room is the preferred location for insertion; however, bedside insertion is acceptable if sterile technique is followed and appropriate diagnostic and therapeutic response to potential complications is accessible.

1. Provide a sterile operative field: use sterile drapes, instruments, and accessories. Perform surgical scrub. Wear gown, cap, gloves, and mask. Have the patient and all personnel in attendance wear a mask.
2. Place the patient in a supine position and expose the upper neck, chest or the groin on the side to be accessed.

- **For subclavian and jugular insertion:** Turn the patient's head slightly to the side to expose the insertion site. The Trendelenberg position may facilitate insertion and prevent air embolism or blood loss.
 - **For femoral insertion:** Flex the patient's knee on the same side as the insertion site. Abduct the thigh on the same side and place the foot across the opposite leg.
3. Shave the access site (optional) and scrub the area with an appropriate antiseptic solution. Isolate the access site with sterile drapes.
 4. Attach a 5 mL, 6 mL, or 10 mL syringe of heparinized saline to each adapter of the catheter. Fill the catheter with 3 to 4 mL heparinized sterile normal saline and clamp immediately. Leave the syringes attached to the adapters.
WARNING: To prevent air embolism, keep the catheter clamped at all times when not attached to a syringe, IV tubing, or bloodlines.
 5. Administer local anesthetic to the skin and underlying tissue at the insertion site.
WARNING: Bending the needles for anesthetic administration is not recommended as it may damage or break the needle.
 6. Flush an 18 G introducer needle with heparinized saline. Insert the needle into the vein in the direction of the blood flow. Aspirate a small amount of blood to ensure the needle is correctly positioned in the vein.
CAUTION: If arterial blood is aspirated, remove the needle and apply immediate pressure to the site for at least 15 minutes. Ensure that arterial bleeding has stopped and that no hematomas have developed before attempting to cannulate the vein again.
 7. a. Disconnect the syringe from the needle and promptly insert the flexible "J" end of the guidewire through the introducer needle. Failure to insert promptly may lead to air embolus or blood loss through the needle. Advance the guidewire into the vein.

CAUTION: For subclavian and jugular insertion, the length of wire inserted is determined by the size of the patient. Cardiac arrhythmia may result if the guidewire passes into the right atrium. If symptoms occur, pull back the guidewire until they are eliminated. If the guidewire meets resistance, do not pull it back through the needle (Figure 1). Remove the guidewire and the needle as a unit, then begin again with new needle and guidewire.


Figure 1. Do not pull guidewire forcibly back through any component.

- b. Withdraw the introducer needle while holding the guidewire firmly in place.
CAUTION: Do not let the guidewire move farther into the vein during the following steps:
8. Make a small incision (0.5 cm) near the guidewire at the skin entry site to facilitate passage of the dilator and catheter.
NOTE: A curved surgical clamp may be used to separate the subcutaneous tissue ("blunt dissection") and create a pocket for the catheter to rest.
 9. Thread a dilator over the end of the guidewire and, with a rotating motion, advance the dilator through the skin and soft tissue until it is just inside the vein. If the dilator does not reduce the resistance sufficiently for the catheter to pass over the wire, progress to the next French size dilator. Once the tract is dilated, remove the tissue dilator and discard.
 10. Thread the tip of the catheter over the guidewire. Open the clamp on the medial extension (white adapter) and remove the syringe to permit the guidewire to exit. Using a rotational motion, advance the catheter through the soft tissue into the vein.
CAUTION: Observe the patient carefully for signs and symptoms of cardiac arrhythmia caused by passage of the

catheter into the right atrium. If symptoms appear, pull back the tip until they are eliminated.

For jugular and subclavian insertions, use X-ray or fluoroscopy to ensure that the tip of the catheter lies at the junction of the superior vena cava and the right atrium.

If using the catheter for hemodialysis, CRRT or apheresis, orient the blue adapter cephalad (Figure 2). This positions the “arterial” inlet away from the wall of the superior vena cava, reducing the potential for one-way inlet obstruction.


Figure 2. Implanted MAHURKAR Elite Catheter (with Curved Extensions)

11. Holding the catheter firmly in place, gently withdraw the guidewire from the lumen. If the guidewire meets more than slight resistance, do not pull it through the catheter. Remove the catheter and the guidewire together as a unit and begin again with new catheter and catheter insertion components.
12. Verify catheter patency by aspirating blood through all three adapters of the catheter. With patency confirmed, inject sterile normal saline into each lumen, followed by appropriate heparinized saline priming volume as indicated on the catheter clamp. Clamp each extension immediately and attach a sterile sealing cap to each adapter.
13. Use X-ray or fluoroscopy to verify correct placement of the catheter tip.
14. Suture the catheter to the skin using the rotatable suture wing.
CAUTION: Do not suture the catheter tubing.
15. Clean the skin around the catheter. Apply dressing to the exit site. Leave the extensions, clamps, adapters, and caps exposed for access by the staff.

SHARP SAFETY DISPOSAL

Discard all procedural components according to hospital protocol in an appropriate biohazard and/or sharps disposal container following applicable regulations.

HEPARINIZATION

To maintain patency between dialysis treatments, keep the lumina of the catheter filled with the appropriate concentration and volume of heparinized saline. In most cases, 5,000 units/ mL is most successful (refer to the priming volume of the catheter). Approved heparinized saline concentrations vary with each institution. Be sure to use those concentrations approved by your facility.

CAUTION: Do not use alcohol or any solution containing alcohol to lock the catheter between dialysis treatments.

Heparinize after use. Before initiating treatment, aspirate indwelling heparinized saline and discard. After treatment, flush well and instill fresh heparinized saline. If the interdialytic period is less than two days, a lower concentration of heparin may be desirable.

In all cases, the patient’s condition must be considered when choosing a heparin regime. Use less heparin in children and in adults with bleeding disorders.

NOTE: Always defer to the physician’s experience and judgment or institution protocol.

Priming Volumes

Catheter		Lumen		
Extensions	Length	Arterial, Red	Venous, Blue	Medial, White
Curved	13 cm	1.2 mL	1.2 mL	0.4 mL
Curved	16 cm	1.3 mL	1.3 mL	0.4 mL
Curved	20 cm	1.4 mL	1.4 mL	0.4 mL
Curved	24 cm	1.5 mL	1.5 mL	0.5 mL
Straight	13 cm	1.2 mL	1.2 mL	0.4 mL
Straight	16 cm	1.3 mL	1.3 mL	0.4 mL
Straight	20 cm	1.4 mL	1.4 mL	0.4 mL
Straight	24 cm	1.5 mL	1.5 mL	0.5 mL
Straight	30 cm	1.7 mL	1.7 mL	0.5 mL

Supplies

- (2) 10–20 mL syringes
- (1) 3 mL syringe
- (3) 20 G (0.9 mm) x 2.5 cm needles
- (1) vial heparin (in concentration approved by your institution)
- (1) vial sterile normal saline
- An appropriate antiseptic.

PREPARATION

1. Prepare supplies on a clean surface.
2. Wash hands thoroughly with soap and water.
3. Scrub the area surrounding the cap and catheter for 5 minutes with an appropriate antiseptic swab. **Allow to air dry.**
4. Open syringe and needle packages. Place the needle on the sterile syringe, using aseptic technique.
5. Remove the tops of the saline and heparin vials and swab the injection area with an appropriate antiseptic. **Allow to air dry.**
6. Prepare the appropriate heparinized saline solution.

PROCEDURE

1. Aspirate indwelling heparinized saline from the catheter before infusing fresh heparinized saline or initiating treatment.
2. Flush each lumen with 10 to 20 mL sterile normal saline.
CAUTION: Before flushing, pull the plunger back to verify blood flow and to ensure there are no blood clots. Do not flush clots through the catheter (see “Thrombus Formation”).
3. Infuse fresh heparinized saline, flushing quickly to ensure that heparinized saline reaches the distal end of the lumen, and clamp immediately. Perform procedure for all three lumens.

Do not infuse against a closed clamp or forcibly infuse a blocked catheter.

Once the lumina have been primed, keep the extensions clamped when not attached to bloodlines or a syringe. If an extension is unclamped, there is high risk of blood loss or air embolism. It also leads to a slightly increased priming volume as a result of the tube returning to its “normal” unclamped state. This creates a vacuum at the tip, causing blood to be drawn into the distal portion of the catheter, potentially resulting in a thrombus.

MANAGEMENT OF ONE-WAY OBSTRUCTION

Suspect one-way obstruction when a lumen can be flushed but blood cannot be aspirated. This is usually caused by tip malposition. Some signs of obstruction:

- There are air bubbles in the tubing set; blood is foamy.
- The venous drip chamber has collapsed or is at a lower level than normal.

One of the following adjustments may eliminate the obstruction:

- Have the patient hold his arms above his head and cough.
- Reposition the patient.
- Flush with saline to push the catheter away from the vessel wall (provided there is no resistance).

- Rotate the catheter to position the arterial inlet away from the wall of the superior vena cava and allow free flow of blood into the arterial lumen.

CAUTION: Do not insert the catheter further into the vein.

- Reverse the blood lines. If the previous methods fail to eliminate the obstruction, connect the arterial bloodline to the venous adapter and the venous bloodline to the arterial adapter. Due to the symmetrical side slot design of the MAHURKAR Elite catheter, similar recirculation is expected in forward and reverse flows.
- Administer thrombolytic agent per hospital protocol and manufacturers' dosing recommendations.

NOTE: Always defer to the physician's experience and judgment or institution protocol.

THROMBUS FORMATION

WARNING: Only a physician or other health care practitioner authorized by and under the direction of such physician should attempt the following procedure.

NEVER FORCIBLY FLUSH AN OBSTRUCTED LUMEN.

If the catheter is blocked, first check to be sure it is not kinked. If thrombus is present, attempt to aspirate the clot gently with a 10 mL syringe.

If infusion is still slow, or blood cannot be withdrawn at all, the physician may choose to dissolve the clot with a thrombolytic agent.

Streptokinase is not recommended; it is reported to be anaphylactogenic in some patients.

The cleared catheter may be used immediately.

NOTE: A catheter that does not clear after a repeated dose of thrombolytic agent may be blocked by a substance other than a blood clot.

PRESSURE INJECTION OF CONTRAST MEDIA PROCEDURE:

1. Maintain sterile technique when handling or using the catheter. Scrub the adapters, sealing caps, clamps, extensions and Y-connector of the catheter with the appropriate antiseptic. Allow a minimum of 5 minutes exposure time to an appropriate antiseptic prior to removing the medial lumen sealing cap.
2. Obtain a visual image to confirm the catheter's tip placement prior to pressure injection of contrast media.
3. Remove the medial lumen sealing cap.
4. Aspirate the indwelling heparinized saline from the medial lumen of the catheter before starting the pressure injection of contrast media.
5. Aspirate the medial lumen for adequate blood return and flush the medial lumen with 10 to 20 mL of normal saline.
6. Warm the contrast media to body temperature (37° C) prior to injection.

WARNING: Failure to warm the contrast media may result in catheter failure.

7. Attach the pressure injection device to the medial lumen per the manufacturer's recommendations.
8. Pressure inject the contrast media according to hospital protocol.

WARNING: Do not exceed the maximum indicated flow rate of 5 mL/sec.

9. Disconnect the pressure injector device from the catheter.
10. Flush the medial lumen with 10 to 20 mL of normal saline then lock the medial lumen per hospital protocol.
11. Immediately clamp the medial lumen and attach a sterile sealing cap.

CENTRAL VENOUS PRESSURE MONITORING (CVP READING):

1. Maintain sterile technique when handling or using the catheter. Scrub the adapters, sealing caps, clamps, extensions and Y-connector of the catheter with the appropriate antiseptic. Allow a minimum of 5 minutes exposure time to an appropriate antiseptic prior to removing the medial lumen sealing cap.
2. Obtain a visual image to confirm the catheter's tip placement prior to central venous pressure monitoring.
3. CVP monitoring is intended to be performed through the medial lumen. Remove the medial lumen sealing cap.

4. Aspirate the indwelling heparinized saline from the medial lumen of the catheter before starting the CVP reading procedure.
5. Ensure that the pressure transducer is at the same level as the right atrium.

CAUTION: A continuous infusion of 3 mL/hr of normal saline is recommended while measuring CVP to improve the accuracy of the CVP reading.

6. Perform the CVP procedure per hospital protocol.

WARNING: CVP readings should not be performed during hemodialysis, apheresis, infusion or pressure injection of contrast media.

7. Flush the medial lumen with 10 to 20 mL normal saline then lock the medial lumen per hospital protocol.
8. Immediately clamp the medial lumen and attach a sterile sealing cap.

"ON/OFF" PROCEDURE

"ON" Procedure

For priming volumes see HEPARINIZATION section of these instructions.

1. Prepare machine and prime bloodlines in the usual manner. Maintain sterile technique when handling or using the catheter. Scrub the adapters, sealing caps, clamps, extensions, and Y-connector of the catheter with an appropriate antiseptic. Allow a minimum of five minutes exposure time to an appropriate antiseptic prior to removing the sealing caps.
2. a. Verify that the "arterial" extension (red adapter) is clamped before removing the sealing cap.
b. Remove the sealing cap from the adapter and attach a luer lock syringe.
c. Confirm patency of the lumen by aspirating indwelling heparinized saline until blood appears. If anticoagulant prime is used, it may be administered at this time.
d. Close the clamp on the extension.
3. Repeat Steps 2a through 2d for the "venous" lumen (blue adapter).

WARNING: When connecting bloodlines to the catheter, do not permit air to enter the blood path. Do not overtighten connections.

4. For hemodialysis or CRRT: Remove the syringe and connect the arterial bloodline to the arterial adapter. Open the clamp on the arterial extension and turn on the blood pump. Prime the extracorporeal circuit completely with the patient's blood and turn off the blood pump. Ensure the venous extension is clamped, then remove the syringe and connect the venous bloodline to the venous adapter. Open the clamp on the venous extension and turn on the pump.

For apheresis: Ensure the extension is clamped, then remove the syringe and connect the arterial bloodline to the arterial adapter; repeat for the venous lumen. Open both clamps and turn on the blood pump.

"OFF" Procedure

1. Observe sterile technique. Stop the blood pump. Close the clamp on the arterial extension and clamp the arterial bloodline at the connection site. Disconnect the arterial bloodline from the adapter of the catheter.
2. Connect a 10-20 mL syringe filled with sterile normal saline to the arterial adapter; open the clamp on the arterial extension and flush the blood from the arterial lumen of the catheter. Reclamp the extension, leaving the syringe attached to the adapter.
3. After the patient's blood has been rinsed back, turn off the blood pump. Clamp the venous extension and disconnect the venous bloodline from the venous adapter of the catheter. Connect a syringe filled with sterile normal saline to the venous adapter. Open the clamp and flush the venous lumen to remove any remaining blood. Reclamp.
4. Connect a syringe containing the appropriate volume and concentration of heparinized saline solution to the venous adapter. Open the clamp on the venous extension and quickly infuse the heparinized saline; reclamp immediately. Remove the syringe and attach a sterile sealing cap to the adapter.
5. Verify that the arterial extension is clamped. Connect a syringe containing the appropriate volume/concentration of heparinized saline solution to the arterial adapter of the catheter. Open the clamp on the arterial extension and

quickly infuse the heparinized saline; reclamp immediately. Remove the syringe and attach a sterile sealing cap to the adapter.

WARNING: Keep the catheter clamped at all times except when connected to bloodlines or syringe during treatment.

EXIT SITE CARE

Keep exit site dry at all times. To reduce the risk of infection while cleaning and dressing the exit site, wash hands thoroughly and wear sterile gloves. Use aseptic technique when handling the catheter and supplies.

Before removing the sealing cap or any connecting line, clamp the catheter and scrub the area surrounding the cap and catheter with antiseptic.

SUPPLIES

mask
gloves
an appropriate antiseptic
sterile, non-adherent gauze sponges
transparent occlusive dressing

PROCEDURE

Maintain aseptic technique when handling or using the catheter. Before beginning, don mask, then wash hands thoroughly. Don non-sterile gloves if desired.

1. Remove the wound dressing carefully from the catheter and exit site.
2. Examine the exit site and surrounding area for inflammation, reddening, or discharge. Using a sterile gauze sponge, palpate the area surrounding the exit site for evidence of tenderness. If signs or symptoms of infection are present, notify the physician immediately. Before cleaning an infected site, verify whether the physician would like the exudates cultured. If so, collect the specimen before continuing.
3. Wash hands again, then don sterile gloves.
4. Clean the exit site in a circular motion from the catheter out, using an appropriate antiseptic.

Allow to air dry. Apply with either sterile swabs or sterile gauze sponges.

CAUTION: Do not use acetone or any solution containing acetone on any part of the catheter. Exposure to these agents may cause catheter damage.

Any ONE of the following can be used as a cleaning agent: aqueous-based povidone-iodine, chlorhexidine, amukin 50%, bacitracin ointment, hydrogen peroxide, Betadine™, Exsept™, 70% isopropyl alcohol, Chloraprep™, Neosporin™, Octenisept™, Dilutis™, Polysporin™, Bactroban™ cream, aqueous chlorhexidine gluconate or chlorhexidine 2% with alcohol. PLEASE NOTE: mixing of these solutions has not been tested and is not recommended.

5. Verify that the sutures are secure in the suture wing. The catheter should not be able to move in and out of the exit site.
6. Apply the transparent occlusive dressing. If desired, sterile gauze sponges can be placed around the catheter at the exit site before the dressing is applied; however, the sponges will impair visual examination of the site. Leave the extensions, clamps, adapters, and sealing caps exposed for access.

CATHETER REPLACEMENT

Per KDOQI guidelines, remove the catheter as soon as it is no longer needed. It is recommended that internal jugular catheters should be used no more than 1 week and femoral catheter should be used for no more than 5 days in bed-bound patients with good exit site care.

Replace the catheter sooner if infection occurs or if there is a progressive increase in venous resistance or progressive decrease in flow rates during the course of hemodialysis, CRRT or apheresis treatment.

NOTE: These recommendations are not intended to supersede the physician's experience and judgment in treating specific patients.

1. Place the patient in a supine position and expose the exit site.
2. Maintain sterile technique. Remove the dressing and examine the exit site and surrounding area for signs or symptoms of infection. If present, the physician must determine whether to postpone catheter replacement until the infection has been successfully treated, or to access the patient in another site.

3. Gown; don sterile gloves; scrub the external portions of the catheter and surrounding area with an appropriate antiseptic, then isolate the site with sterile drapes.
4. Flush all three lumina of the new catheter with 3 to 4 mL heparinized saline solution and clamp immediately.
5. Cut and remove the old sutures carefully from the skin.
6. Ensure that the clamp on the medial extension (white adapter) of the indwelling catheter is closed. Simultaneously remove the sealing cap from the white medial adapter while placing the "J" end of the guidewire into the medial lumen of the adapter. Open the clamp and pass the guidewire through the catheter to the appropriate position.
CAUTION: For subclavian and jugular insertion, observe the patient for cardiac arrhythmia which may result if the guidewire passes into the right atrium.
7. Hold the guidewire steady to keep it from exiting the vein. Carefully remove the catheter by sliding it over the wire. Discard the old catheter.
8. Pass the new catheter over the guidewire following steps 10 through 15 in INSERTING THE CATHETER.

REFERENCES

Seldinger, S.I. "Catheter Replacement of Needle in Percutaneous Arteriography: New Technique." ACTA Radiologica 39 (1953): 368-76.
Coe, A.J. et al: "Triple lumen catheter" British Journal of Hospital Medicine, April 1988: 313-316.
KDOQI guidelines, 2006

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Caution, consult
accompanying
documents

Non-pyrogenic	Pyrogenfri
Apyrogene	Μη πυρετογόνα
Pyrogenfrei	Αpyρογεννή
Apirogene	Nem pirogén
Apirógeno	Свободные от
Pyrogenfri	пирогенов
Niet-pyrogene	Produkt
Apirogénico	niegorączkotwórczy
Pyrogeniton	Pirojenik değildir

Keep away
from sunlight

0123

IFU No. PT00047070

MAHURKAR is a trademark of Sakharam D. Mahurkar, MD, used under license.

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG.

™*Trademark of its respective owner.

©2012 Covidien.

Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.

EC REP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.

www.covidien.com

REV 05/2016

COVIDIEN 15 Hampshire Street, Mansfield, MA 02048

Job #: <small>(If applicable)</small>	VT716	Business Unit:	Patient Monitoring & Recovery/ Dialysis	Colors:
Revision Date:	05/06/2016	Artwork/Part #:	PT00047070	 Black
File Name:	PT00047070_GNRC_IFU			<input type="checkbox"/> PMS XXX
Template ID:	N/A			<input type="checkbox"/> PMS XXX
Dimensions:	3.75" W x 11.0" H	Dieline Name: <small>(If applicable)</small>	N/A	<input type="checkbox"/> PMS XXX
Print Process:	Offset	Package/Item Type:	IFU	<input type="checkbox"/> PMS XXX
Designer(s):	N. Gouveia	Application/Version:	InDesign CC / 2015	<input type="checkbox"/> PMS XXX
Associated to: <small>(If applicable)</small>	SEE ATTACHMENTS TAB IN AGILE			
Comments/ Printer Notes:	N/A			