# Palindrome™

## Precision Chronic Catheter Symmetrical Tip

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Sterilized with ethylene oxide. Sterile and non-pyrogenic in unopened and undamaged package.

**CAUTION:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

READ ALL INSTRUCTIONS, WARNINGS,  
AND CAUTIONS CAREFULLY PRIOR TO USE.

DESCRIPTION:

The Palindrome™ precision 14.5 Fr/Ch chronic catheter is a radiopaque, urethane catheter with felt cuff and dual extensions.

Each extension has an in-line clamp and luer-lock adapter that is color coded: red for “arterial” outflow of blood, blue for “venous” return. Sealing caps are included with the catheter.

Dimensions:

CATHETER O. D.	14.5 Fr/Ch (4.8 mm)					
OVERALL LENGTH	36 cm	40 cm	45 cm	50 cm	61 cm	72 cm
IMPLANT LENGTH	19 cm	23 cm	28 cm	33 cm	44 cm	55 cm

PALINDROME™ PRECISION CATHETER FLOW TABLE  
WITH SIDE SLOTS:

FLOW RATE (mL/min)	Implant Lengths											
	19 cm		23 cm		28 cm		33 cm		44 cm		55 cm	
	+VP*	-AP**	+VP*	-AP**	+VP*	-AP**	+VP*	-AP**	+VP*	-AP**	+VP*	-AP**
200	18	-19	21	-23	23	-24	26	-28	33	-34	40	-41
250	24	-28	29	-32	33	-35	37	-40	45	-49	54	-58
300	33	-38	39	-44	45	-48	51	-54	58	-66	70	-78
350	44	-50	50	-58	59	-64	67	-72	74	-86	88	-102
400	56	-64	65	-74	75	-82	83	-93	93	-109	111	-130
450	71	-80	79	-92	91	-103	102	-117	112	-136	135	-164
500	87	-98	96	-112	108	-126	120	-143	133	-167	161	-203

\*postive (+VP) Venous Pressure (mm/Hg)  
\*\*negative (-AP) Arterial Pressure (mm/Hg)

Recirculation rates were tested on a bench top with a flow rate of 400 mL/min. Results were 1.5% and 3.1% in forward and reverse directions respectively.

Dynamic flow rates were measured with deionized water.

MR SAFE

MRI testing demonstrated that the materials used to make the catheter were non-magnetic, non-metallic, and non-conducting. Therefore, based on the findings from the MRI tests and the materials used to make this device, it is considered to be MR-Safe according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

INDICATIONS:

The Palindrome™ precision 14.5 Fr/Ch chronic catheter is intended for acute and chronic hemodialysis, apheresis, and infusion. It may be inserted either percutaneously or by cutdown. Catheters greater than 40 cm implant length are indicated for femoral insertion.

CONTRAINDICATIONS:

Do not use this catheter in thrombosed vessels or for subclavian puncture when ventilator is in use.

POTENTIAL COMPLICATIONS:

- sepsis  
exit site infection  
air embolism  
hemorrhage  
pneumothorax  
cardiac tamponade  
pulmonary emboli  
catheter thrombosis  
femoral nerve damage  
femoral artery dissection  
lower extremity ischemia  
arterial puncture  
deep vein thrombosis of the lower extremity
- thrombosis/stenosis of vein  
cardiac arrhythmia  
subcutaneous tunnel infection  
hemothorax  
hematoma  
trauma to major vessel or right atrium  
brachial plexus injury  
retroperitoneal bleed  
femoral artery damage  
femoral vein occlusion  
pulmonary embolism  
mediastinal widening  
endocarditis

## WARNINGS AND PRECAUTIONS:

- The catheter should be inserted and removed only by a qualified, licensed physician or other healthcare practitioner authorized by and under the direction of such physician.
- The medical techniques and procedures described in these instructions do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.
- Observe sterile technique at all times when handling or using the catheter.
- Do not use the catheter if package has been previously opened or is damaged.
- Do not use the catheter if it appears damaged or defective.
- Avoid air embolism by keeping catheter extension tubing clamped at all times when not in use and by filling the catheter with sterile saline prior to implantation. With each tubing change, purge air from the tubing and aspirate any air in catheter.
- Do not clamp extension tubing with insertion stylets loaded as it may result in stylet damage.
- Percutaneous insertion of the catheter in the subclavian vein can be technically difficult. The right internal jugular is preferable.
- To avoid vessel perforation and damage, do not insert the guidewire, dilators, or valved pull-apart sheath/introducer forcibly.
- Do not insert the valved pull-apart sheath/introducer further than necessary: depending upon patient size and access site, it may not be necessary to insert the entire length of the introducer into the vessel.
- The valved pull-apart sheath/introducer is designed to reduce blood loss and the risk of air intake but it is not a hemostasis valve.
- The valved pull-apart sheath/introducer is not intended to create a complete two-way seal nor is it intended for arterial use.
  - The valve will substantially reduce air intake. At -12 mm Hg vacuum pressure the valved pull-apart sheath/introducer may allow up to 4 mL/sec of air to pass through the valve.
  - The valve will substantially reduce the rate of blood flow but some blood loss through the valve may occur.
- Use the guidewire straightener to insert the "J" end of the guidewire into the introducer needle. Do not insert or withdraw the guidewire forcibly from any component; the wire could break or unravel.
- Do not nick the catheter when suturing.
- Do not tie the suture too tightly at the venotomy site.
- Prolonged exposure to ultraviolet light can damage the catheter.
- Do not use acetone on any part of the catheter. Aqueous-based povidone iodine, ExSept™, Hibiclens™ (Chlorohexidine), amukin 50%, hydrogen peroxide, Neosporin™ antibiotic ointment, bacitracin ointment, Bactroban™ cream, isopropyl alcohol 70%, ChloroPrep™ can be used. Inter-mixing of these solutions has not been tested and is not recommended.
- Overtightening catheter connections can crack some adapters.
- Do not clamp the dual lumen portion of the catheter; clamp only the extensions. Use only smooth-jawed forceps for clamping when not using the clamp supplied with the catheter.
- Clamping the catheter repeatedly in the same spot could weaken the tubing: change the position of the clamp regularly to prolong the life of the tubing. Avoid clamping near the adapter and hub.
- Exercise caution when using sharp instruments near the catheter. Catheter tubing can tear when subjected to nicks, excessive force, or rough edges.
- Inspect the catheter frequently for nicks, scrapes, cuts, etc. which could impair its performance.
- When infusing heparin, flush quickly and clamp immediately to ensure heparin reaches the distal end of the lumen. Do not infuse against a closed clamp or forcibly infuse a blocked catheter: back pressure could force the adapter out of the tubing.
- Remove the catheter as soon as it is no longer necessary.
- When removing the catheter, DO NOT use a sharp, jerking motion or undue force; this may tear the catheter. Free the cuff and surfaces from the tissue prior to removal.
- Discard the catheter after single use. Do not resterilize.
- This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

## RECOMMENDED

- Use only luer-lock (threaded) connections (including syringes, bloodlines, and sealing caps) with the catheter's adapters.
- Use a straight tip guidewire when removing a Palindrome™ precision catheter with side slots for an over-the-guidewire catheter exchange.

INSERTING THE PALINDROME™ PRECISION 14.5 FR/CH  
CHRONIC CATHETER

Sterile Supplies Required

catheter	normal saline
prepping agents	#11 scalpel
drapes	tunneler
mask, gloves, gown	needle holder
syringes/needles	suture with curved needle
local anesthetic	sealing caps
sponges	wound dressing
heparin (in concentrations approved by your institution)	razor (optional)

Percutaneous placement also requires the following items:

- 18 G introducer needle
- 0.038" (0.97 mm) J/straight stainless steel guidewire
- 16 Fr/Ch pull-apart sheath/introducer
- 12 mL syringe
- tissue dilators (optional)
- insertion stylets (optional)

INSERTION SITE

The Palindrome™ precision 14.5 Fr/Ch chronic catheter is ideally placed in the right atrium via the right internal jugular. While the catheter also can be placed in the external jugular, subclavian, femoral, or saphenous vein, the right internal jugular is strongly recommended for the following reasons:

- The internal jugular permits easier positioning of the catheter tip in the right atrium.
- The size and location of the external jugular makes insertion difficult.
- Use of the subclavian vein for catheter placement may result in subclavian vein stenosis. Subclavian vein stenosis may prevent the future use of the ipsilateral extremity for permanent access.
- Subclavian vein placements are at risk for higher insertion complications.
- The saphenous vein should be used only as a last alternative due to the possibility of insertion complications.
- The risk of thrombosis, infection and bleeding are increased with femoral placements.

For optimal catheter function: Per Dialysis Outcome Quality Initiatives (DOQI), the catheter tip must be adjusted to the level of the caval atrial junction or beyond to ensure optimal blood flow. For this reason, it is preferable to insert the catheter on the patient’s right side. For large patients, and for patients whose right-side veins are unusable, 40 cm, 45 cm, 50 cm catheters are available for placement via the left jugular or subclavian vein. The 61 cm and 72 cm length is intended for femoral placement. The catheter tip in femoral placement should be above the common iliac bifurcation, within the inferior vena cava.

PREPARATION

The operating room or IR suite are the preferred locations for catheter placement. **Both cutdown and percutaneous procedures require confirmation of correct placement by fluoroscopy or chest x-ray.**

1. Provide a sterile operative field: use sterile drapes, instruments, and accessories. Perform surgical scrub. Wear a gown, cap, gloves, and mask. Have the patient wear a mask.
2. Place the patient in a supine position and expose the upper chest or the groin side to be accessed.

**For subclavian and jugular insertion:** Turn the patient’s head slightly to the side to expose the insertion site. The Trendelenberg position may facilitate insertion.

**For femoral insertion:** Flex the patient’s knee on the same side as the insertion site. Abduct the thigh on the same side and place the foot across the opposite leg.

3. Shave the access site (optional) and prep the area in the established manner. Isolate the access site with sterile drapes.

NOTE: When performing a subclavian insertion on a patient with large breasts, it is best to draw landmarks while the patient is sitting up to prevent catheter tip migration.

4. Fill the catheter with sterile heparinized saline and clamp the extensions immediately.
5. Administer local anesthetic to the skin and underlying tissue at the insertion site.

**WARNING: To prevent air embolism, keep the catheter clamped at all times when not attached to a syringe, IV tubing, or bloodlines.**

**NOTE: If using insertion stylets, catheter will not be filled until after placement. Load stylets into catheter lumens such that a minimum length of 8 cm extends distal from the catheter adapters, in order to allow for proper catheter tunneling.**

## PERCUTANEOUS PROCEDURE MODIFIED SELDINGER (LITTLEFORD-SPECTOR) TECHNIQUE

### CANNULATING THE VESSEL

1. Flush an 18 G introducer needle with heparinized normal saline. Insert the needle through the primary incision and advance it into the vein, in the direction of blood flow. Aspirate a small amount of blood to ensure the needle is correctly positioned in the vein.


**CAUTION: If arterial blood is aspirated, remove the needle and apply immediate pressure to the site for at least 15 minutes. Ensure that the bleeding has stopped and that no hematomas have developed before attempting to cannulate the vein again.**

- 2a. Disconnect the syringe from the needle and promptly insert the flexible "J" end of the guidewire through the introducer needle. Failure to insert the wire promptly may lead to blood loss through the needle. Advance the wire into the vein.

**CAUTION FOR JUGULAR AND SUBCLAVIAN INSERTION: The length of wire inserted is determined by the size of the patient. Cardiac arrhythmia may result if the guidewire passes into the right atrium. If symptoms occur, pull back the guidewire until they disappear.**

**If the guidewire meets resistance, do not pull it back through the needle. Remove the wire and the needle together as a unit and begin again with new needle and guidewire (Figure 1).**

NOTE: If utilizing the stylets for insertion, the guidewire provided is recommended. Otherwise, greater than a 0.035" hydrophilic or greater than a 0.038" stainless steel wire is contraindicated.


- 2b. Withdraw the introducer needle, leaving the guidewire in the vein.

### CREATING THE SUBCUTANEOUS TUNNEL

NOTE: For ease in dressing the exit site and for patient comfort, locate the subcutaneous tunnel below the entry site. A tunnel with a wide, gentle arc lessens the risk of kinking the catheter at the cuff. The tunnel should be short enough to keep the Y-hub of the catheter approximately 3 cm outside the exit site, yet long enough to keep the cuff a minimum of 2 cm deep, and a minimum of 3 cm from the insertion site. For femoral placement it is recommended to tunnel laterally and superiorly from the access site toward the pelvic region, away from the groin. This reduces the likelihood of infections and allows easy access to the catheter for dialysis.

1. Make a small incision at the insertion site. (The right internal jugular is the preferable site for percutaneous placement as subclavian placement is technically difficult.) Make a second incision parallel to the first at the exit site. Make the exit incision just long enough to accommodate the cuff, approximately 1 cm.
2. Use blunt dissection to create the subcutaneous tunnel.
  - a. Attach the catheter to the tunneler by sliding the catheter tip onto the bifurcated tunneler tines until the catheter tip meets the base of the tines. Slide the sheath completely over the connection until it stops, being careful that the sheath smoothly transitions over the catheter tip (see Figure 2). If desired, bend the tunneler into a wide arc to make a curved tunnel.

- b. Insert the tunneler into the exit site and create a short subcutaneous tunnel, emerging at the insertion site. The catheter will thread through the tissue as the tunnel is created. Ensure the catheter passes through the tunnel to the primary insertion site.
- c. Remove the catheter carefully from the tines. To do so, slide the sheath back, grasp the tip of the catheter, and gently pull the catheter from the tunneler tines. **Discard the tunneler.**

## INSERTING THE CATHETER USING THE VALVED PULL-APART SHEATH/INTRODUCER

OPTIONAL: To ease insertion of the valved pull-apart sheath/ introducer, the vein can be pre-dilated with the dilator(s) provided.

1. Thread the dilator(s) over the end of the guidewire and advance it into the vein using a rotating motion to assist passage through the tissue.

**CAUTION: Do not force the dilator(s). Ensure that the guidewire does not advance further into the vein.**

Remove and discard the dilator(s).

2. Remove the valved pull-apart sheath/introducer assembly from the package; remove dilator from the sheath assembly and discard the clear sheath guard. (Figure 3).

Figure 3

3. Insert the dilator through the valve sheath opening and lock in place using the rotating red hub. (Refer to Figures 4, not locked position and 5 the locked position.)


Figure 4


Figure 5

4. a. Thread the locked valved pull-apart sheath/introducer assembly over the end of the guidewire.

**CAUTION: To avoid damaging the tissue and the sheath tip, do not let the sheath advance over the dilator. The two must be grasped as one unit.**

- b. With a rotational motion, advance the assembly into the vein only as far as necessary. Do not force the introducer into the vessel. Do not insert it further than necessary for the patient's size and access site. Ensure that the guidewire does not move further into the vein.
  - c. Hold the sheath in place, unlock the rotating red hub, and gently remove the dilator and guidewire together (discard the dilator and the guidewire).
5. Advance catheter through the valve. To prevent kinking the catheter, it may be necessary to advance in small steps by grasping the catheter close to the sheath (Figure 6).

**WARNING: Ensure the catheter is filled with heparinized saline and is free of air bubbles before inserting it into the vein.**



Figure 6

6. Continue inserting the catheter through the valved sheath and into the vein. For optimal catheter function: Per Kidney Dialysis Outcome Quality Initiatives (K/DOQI), the catheter tip must be adjusted to the level of the caval atrial junction or beyond to ensure optimal blood flow (Figure 8).
7. Aspirate to verify patency and clamp the extension.


8. After the catheter has been positioned, grasp both tabs of the sheath firmly and pull apart to crack the sheath handle in half.  
NOTE: Steps 9, 10 and 11 refer to Figure 7.
9. Once both tabs are split, note that the valve in the sheath has been broken.
10. Peel the sheath partially away from the catheter.  
NOTE: Some resistance will be experienced while pulling the sheath off of the catheter.

11. Remove the sheath from the patient by holding the catheter in place and simultaneously pulling the separated tabs away from the entry site at a 180° angle. The sheath will separate from the catheter, leaving the catheter in the vein.

**CAUTION: Do not allow the catheter to move out of the vein with the sheath. Ensure that the vein is not bleeding around the catheter.**

12. Use fluoroscopy or portable x-ray to view the catheter. Position the catheter correctly in the right atrium.  
Note: the catheter tip in femoral placement should be above the common iliac bifurcation, within the inferior vena cava.

**Figure 8: Implanted Catheter**

To aid in catheter tip visualization a radiopaque marker band has been placed adjacent to the most proximal side slot.


13. Confirm correct placement and catheter function by aspirating venous blood from both the arterial and venous lumina, then flush 5 mL of sterile normal saline into each lumen. Follow with heparinized saline. Clamp the extensions immediately utilizing the positive pressure technique. (See priming volumes under HEPARINIZATION.) Attach a sterile sealing cap to each adapter.  
**WARNING: To prevent air embolism, expel all air from syringes before injecting solutions.**
14. Suture the entry site. The exit site should not require suture.
15. Suture the catheter hub to the skin utilizing a polypropylene suture or equivalent 3-0 or 4-0 suture on the extension wings.  
**CAUTION: Do not suture through any part of the catheter.** Remove skin suture by the 4th or 5th day to diminish erythema in the area. The catheter can also be immobilized with either gauze or transparent dressings.
16. Apply dressing to catheter exit site and insertion site.

## SHEATHLESS INSERTION OR CATHETER EXCHANGE USING VENATRAC™ OVER-THE-WIRE INSERTION STYLETS (NOT APPLICABLE TO FEMORAL PLACEMENT)

1. Load stylets fully into catheter and lock stylet hubs to lumen adapters. Ensure that the white flat side of the longer stylet containing the stylet sidehole is oriented against the septum of the catheter lumen (Figure 9).

NOTE: Flush the VenaTrac™ insertion stylets with sterile saline prior to use. Wetting the outer surfaces of the stylets as well may facilitate advancement into the catheter.

**Figure 9**

**OPTIONAL: To ease insertion of the catheter, the vein can be pre-dilated with the dilator(s) provided.**

2. Thread the dilator(s) over the end of the guidewire and advance it into the vein using a rotating motion to assist passage through the tissue.


**CAUTION: Do not force the dilator(s). Ensure that the guidewire moves freely within the dilator(s) and does not advance further into the vein.**

Remove and discard the dilator(s).

3. Thread the guidewire 1 cm into the distal tip of the longer stylet such that it exits the stylet sidehole (Figure 9). To assist the guidewire exit, the longer stylet tip should be flexed 45°.
4. Continue threading the wire into the distal tip of the adjacent shorter stylet (Figure 9) and advance until the wire passes out the proximal end of the catheter-stylet assembly.
5. While maintaining control of the guidewire, pinch the catheter-stylet slightly, approximately 2 cm from the distal tip, and advance into the vein (Figure 10).

NOTE: For difficult placements, if so desired, the catheter may be rotated in a clockwise direction to ease advancement. If excessive resistance is felt, stop advancement and retract the entire assembly. Verify that the guidewire has not kinked and re-attempt catheter advancement.

**CAUTION: Do not force the assembly. Ensure that the guidewire does not advance further into the vein.**

**Figure 10**

6. Continue to pinch slightly and advance the catheter-stylet assembly until catheter tip is positioned as desired (Figure 10). Per Dialysis Outcome Quality Initiatives (DOQI), the catheter tip must be adjusted to the level of the caval atrial junction or beyond to ensure optimal blood flow (Figure 8).
7. Holding the catheter in place, gently remove the guidewire and discard. Next, immediately unlock, remove and discard each of the two VenaTrac™ insertion stylets.

**CAUTION: Do not remove insertion stylets while guidewire is in place as it may result in catheter damage. Ensure that the catheter does not move out of the vein while removing insertion stylets.**

**CAUTION: If guidewire meets resistance, do not pull it back through the catheter-stylet assembly. Remove the catheter, insertion stylets, and wire together as a unit and begin again with a new catheter.**

8. Aspirate to verify patency and clamp the extension.
9. Use fluoroscopy or portable x-ray to view the catheter. Position the catheter correctly in the right atrium (Figure 8).
10. Confirm correct placement and catheter function by aspirating venous blood from both the arterial and venous lumina, then flush 5 mL of sterile normal saline into each lumen. Follow with heparinized saline. Clamp the extensions immediately utilizing the positive pressure

technique. (See priming volumes under HEPARINIZATION.) Attach a sterile sealing cap to each adapter.

**WARNING: To prevent air embolism, expel all air from syringes before injecting solutions.**

11. Suture the entry site. The exit site should not require suture.
12. Suture the catheter hub to the skin utilizing a polypropylene suture or equivalent 3-0 or 4-0 suture on the extension wings.

**CAUTION: Do not suture through any part of the catheter.** Remove skin suture by the 4th or 5th day to diminish erythema in the area. The catheter can also be immobilized with either gauze or transparent dressings.

13. Apply dressing to catheter exit site and insertion site.

## CUTDOWN PROCEDURE

1. Make an adequate incision following the skin lines over the desired vessel.
2. Ensure the chosen vein is large enough to accept the catheter.
3. If using a large vein such as the jugular, immobilize the vein and place a purse string suture around the venotomy site. If using a smaller vein such as the saphenous, tie down the distal ligature, but leave the proximal one untied until the catheter is in position.
4. Make a counter incision for the exit site of the subcutaneous tunnel. Make the incision just long enough to accommodate the cuff, approximately 1.0 cm.

NOTE: For ease in dressing the exit site and for patient comfort, locate the subcutaneous tunnel below the entry site. A tunnel with a wide, gentle arc lessens the risk of kinking the catheter at the cuff. The tunnel should be short enough to keep the Y-hub of the catheter approximately 3 cm outside the exit site, yet long enough to keep the cuff a minimum of 2 cm deep, and a minimum of 3 cm from the insertion site.

5. Use blunt dissection to create a short subcutaneous tunnel, emerging at the cutdown incision.
  - a. Attach the catheter to the tunneler by sliding the catheter tip into the bifurcated tunneler tines until the catheter tip meets the base of the tines. Slide the sheath completely over the connection until it stops, being careful that the sheath smoothly transitions over the catheter tip (see Figure 2 in Percutaneous Procedure section). If desired, bend the tunneler into a wide arc to make a curved tunnel.
  - b. Insert the tunneler into the exit site and create a short subcutaneous tunnel, emerging at the cutdown incision. The catheter will thread through the tissue as the tunnel is created. Ensure the catheter passes through the tunnel to the cutdown incision.
  - c. Remove the catheter carefully from the tines. To do so, slide the sheath back, grasp the tip of the catheter, and gently pull the catheter from the tunneler tines. **Do not tug the catheter tubing.** Discard the tunneler.
6. Use forceps to grasp the designated venotomy vein transversely across its entire width.
7. Insert a #11 blade at the midpoint of the vein diameter; cut anteriorly for a controlled venotomy. When the vein is released from the forceps, it will display an oval venotomy. This may be dilated further using a mosquito hemostat.

**WARNING: Ensure the catheter is filled with heparinized saline and is free of air bubbles before inserting it into the vein.**

8. Grasp the end of the catheter with forceps and insert it into the vein. Slacken the proximal ligature to allow the catheter to pass with minimal back bleeding. Advance the catheter into the vein: for subclavian and jugular insertion, advance the catheter tip into the right atrium. For optimal catheter function: Per Dialysis Outcome Quality Initiatives (DOQI), the catheter tip must be adjusted to the level of the caval atrial junction or beyond to ensure optimal blood flow. (See Figure 8 in Percutaneous Procedure section.)
9. Pull the purse string suture (or proximal ligature) closed, but do not tie it before determining the exact position of the catheter.
10. Use fluoroscopy or portable x-ray to view the catheter. Move or reposition the catheter until it is positioned correctly in the right atrium as indicated in Figure 8 in Percutaneous Procedure section. Note: the catheter tip in femoral placement should be above the common iliac bifurcation, within the inferior vena cava.
11. Tie the purse string suture (or proximal ligature) snugly around the catheter. Tie it just tightly enough to control bleeding at the venotomy; do not occlude the catheter.

**CAUTION: Do not allow the cuff to enter the vein or the venotomy. Ideally, it should not touch the vein.**

12. Confirm correct placement and catheter function by aspirating venous blood from both the arterial and the venous lumina, then flush 5 mL sterile normal saline into each lumen. Follow with heparinized saline.

Clamp the extensions immediately utilizing the positive pressure technique. (See priming volumes under HEPARINIZATION.) Attach a sterile sealing cap to each adapter.

**WARNING: To prevent air embolism, expel all air from syringes before injecting solutions.**

- 13. Suture the entry site. The exit site should not require suture.
- 14. If a skin suture is desired, suture the catheter hub to the skin utilizing a polypropylene suture or equivalent 3-0 or 4-0 suture on the extension wings.

**CAUTION: Do not suture through any part of the catheter.** Remove skin suture by the 4th or 5th day to diminish erythema in the area. The catheter also can be immobilized with either gauze or transparent dressings.

- 15. Apply dressing to the catheter exit site and cutdown incision.

**REMOVING THE PALINDROME™ PRECISION 14.5 FR/CH CHRONIC CATHETER**

To remove the catheter, free the cuff from the tissue and pull the catheter gently and smoothly. Do not use sharp, jerking motions or undue force. This could tear the catheter.

**WARNING: If the catheter offers resistance, do not pull further. Perform a cutdown and remove all sutures at the venotomy site.**

**HEPARINIZATION**

To maintain patency between dialysis or apheresis treatments, keep the lumina of the catheter filled with the appropriate concentration and volume of heparin. In most cases, 5,000 units/mL is most successful (refer to the priming volume of the catheter). Approved heparin concentrations vary with each institution. Be sure to use those concentrations approved by your facility.

**PRIMING VOLUMES**

Catheter	Overall Length	Arterial	Venous
ADULT	36 cm	1.6 mL	1.6 mL
	40 cm	1.9 mL	1.9 mL
	45 cm	2.1 mL	2.1 mL
	50 cm	2.3 mL	2.3 mL
	61 cm	2.6 mL	2.6 mL
	72 cm	3.1 mL	3.1 mL

Heparinize only after use. Before initiating treatment, aspirate indwelling heparin and discard. After treatment, flush well and instill fresh heparin. If the interdialytic period is less than two days, or if apheresis is performed daily, a lower concentration of heparin may be desirable.

**In all cases, the patient’s condition must be considered when choosing a heparin regime. Use less heparin in children and in adults with bleeding disorders.**

**SUPPLIES**

- 10-20 mL syringes
- 3 mL syringe
- 20 G 2.5 cm needles
- vial heparin (in concentration approved by your institution)
- vial sterile normal saline
- povidone-iodine swabs

**PREPARATION**

- 1. Prepare supplies on a clean surface.
- 2. Wash hands thoroughly with soap and water.
- 3. Scrub the area surrounding the cap and catheter for 5 minutes with a povidone-iodine swab. **Allow to air dry.**
- 4. Open syringe and needle packages. Place the needle on the sterile syringe, using aseptic technique.
- 5. Remove the tops of the saline and heparin vials and swab the injection area with povidone-iodine. **Allow to air dry.**
- 6. Prepare the appropriate dilute heparin solution.

**PROCEDURE**

- 1. Remove the sealing cap and aspirate indwelling heparin from the catheter before infusing fresh heparin or initiating treatment.
- 2. Flush the lumen with 10 to 20 mL sterile normal saline.  
**CAUTION: Before flushing, pull the plunger back to verify blood flow and to ensure there are no blood clots. Do not flush clots through the catheter (see “Thrombus Formation”).**
- 3. Infuse fresh heparin, flushing quickly to ensure that heparin reaches the distal end of the lumen, and clamp immediately. Infusing or clamping too slowly may cause heparin to exit the catheter from the catheter slot, leaving the distal catheter tip unprotected from

thrombus formation. Do not infuse against a closed clamp or forcibly infuse a blocked catheter: back pressure could force the adapter to loosen and potentially come out of the tubing. Perform procedure for both lumina.

Once the lumen has been primed, keep the extension clamped when not attached to a bloodline or syringe. If the extension is unclamped, the priming volume will increase slightly as a result of the tube returning to its “normal” unclamped state. This creates a vacuum at the tip, causing blood to be drawn into the distal portion of the catheter, ultimately resulting in a thrombus.

## MANAGEMENT OF ONE-WAY OBSTRUCTION

One-way obstruction, which exists when a lumen can be flushed easily, but blood cannot be aspirated, usually is caused by tip malposition. One of the following adjustments may resolve the obstruction:

- Reposition the patient.
- Have the patient cough.
- Provided there is no resistance, flush the catheter vigorously with sterile normal saline to try to move the tip away from the vessel wall.

The following procedures may require a physician’s order:

- Consider using a thrombolytic agent.
- If the one-way obstruction exists in the arterial lumen, consider reversing the bloodlines. The patient may be dialyzed by connecting the arterial bloodline to the venous adapter and the venous bloodline to the arterial adapter.

## THROMBUS FORMATION

NEVER FORCIBLY FLUSH AN OBSTRUCTED LUMEN.

If either lumen develops a thrombus, first attempt to aspirate the clot with a syringe. If aspiration fails, the physician may attempt to lyse the clot with a thrombolytic agent.

**WARNING: Thrombolytic agents may cause systemic fibrinolysis if infused into the circulation. Refer to manufacturer’s instructions, indications for use, and contraindications before using thrombolytic agents. Streptokinase is not recommended; it is reported to be anaphylactogenic.**

## DIALYSIS

For priming volumes, see Heparinization section of these instructions.

1. Prepare dialysis machine and prime bloodlines in the usual manner. Maintain sterile technique when handling or using the catheter. Scrub the adapters, sealing caps, clamps, extension tubes, and Y-connector of the catheter with an aqueous-based povidone-iodine solution.
2.
  - a. Verify that the arterial extension is clamped before removing the sealing cap.
  - b. Remove the sealing cap from the adapter and attach a luer-lock syringe.
  - c. Confirm patency of the lumen by aspirating indwelling heparin until venous blood appears. Leave the syringe in place.
  - d. Close the clamp on the extension.
3. Repeat Steps 2a. through 2d. for the venous lumen.

**WARNING: When connecting bloodlines to the catheter, do not permit air to enter the blood path.**

4. Remove the syringe and connect the arterial bloodline to the arterial (red) adapter. Open the clamps on the arterial extension and arterial and venous bloodlines and turn on the blood pump.
5. Prime the extracorporeal circuit with the patient’s blood and turn off the blood pump. Ensure the venous extension is clamped, then remove the syringe and connect the venous bloodline to the catheter’s venous (blue) adapter. Open the clamps on the venous extension and bloodline and turn on the pump.
6. Begin treatment.

## POST DIALYSIS

Prepare syringes with sterile normal saline and heparin.

1. Stop the blood pump. Close the clamp on the arterial extension and clamp the arterial bloodline at the connection site. Disconnect the arterial bloodline from the adapter of the catheter.
2. Connect a 10-20 mL syringe filled with sterile normal saline to the arterial adapter; open the clamp on the arterial extension and flush the blood from the arterial lumen of the catheter. Reclamp the extension, Heparinize the lumen with the appropriate volume/concentration of heparin.
3. Rinse back the blood in the extracorporeal circuit via the catheter’s venous lumen.
4. After the patient’s blood has been rinsed back, turn off the blood pump. Clamp the venous extension and disconnect the venous bloodline from the venous adapter of the catheter.

5. Connect a 10-20 mL syringe filled with sterile normal saline to the venous adapter. Open the clamp on the venous extension and flush all remaining blood from the venous lumen of the catheter. Reclamp. Heparinize the lumen with the appropriate volume/concentration of heparin.
6. Ensure the clamps are closed on both extensions. Remove syringes and attach a sealing cap to each adapter.

**WARNING: Keep the catheter clamped at all times except when connected to bloodlines or syringe during treatment.**

### **CATHETER CARE GUIDES**

For further information and a copy of Catheter Care guidelines for the clinician, contact your Covidien representative. In the United States call 1-800-962-9888 and for international inquiries, call 508-261-8000.

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Caution, consult  
accompanying  
documents

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
**CE**  
0123

IFU No. PT00006783

COVIDIEN, COVIDIEN with logo and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. <sup>TM</sup>\* brands are trademarks of their respective owner.

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REV 04/2019

COVIDIEN 15 Hampshire Street, Mansfield, MA 02048

<b>Job #:</b> <small>(If applicable)</small>	VT764	<b>Business Unit:</b>	RCS	<b>Colors:</b>
<b>Revision Date:</b>	04/12/2019	<b>Artwork/Part #:</b>	PT00006783	<input checked="" type="checkbox"/> BLACK
<b>File Name:</b>	PT00006783_GNRC_IFU			<input type="checkbox"/> PMS XXX
<b>Template ID:</b>	PT00006783_GNRC_IFU			<input type="checkbox"/> PMS XXX
<b>Dimensions:</b>	3.75" W x 11.00" H	<b>Dieline Name:</b> <small>(If applicable)</small>	N/A	<input type="checkbox"/> PMS XXX
<b>Print Process:</b>	Offset	<b>Package/Item Type:</b>	IFU (Instructions)	<input type="checkbox"/> PMS XXX
<b>Designer(s):</b>	J. Frumerie	<b>Application/Version:</b>	Indesign CC (2019)	<input type="checkbox"/> PMS XXX
<b>Associated to:</b> <small>(If applicable)</small>	8888119360P, 8888123400P, 8888128450P, 8888133500P, 8888119364P, 8888123404P, 8888128454P, 8888133504P, 8888145014P, 8888145015P, 8888145016P, 8888145017P, 8888145018P, 8888145039P, 8888145040P, 8888145041P, 8888145042P, 8888146044P			
<b>Comments/ Printer Notes:</b>	CE Symbol Master 8.9			

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© 2012 Covidien.
Covidien Ireland Limited, IDA Business & Technology Park, Tullamore, Ireland.