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
4a BD

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SPECIFICATION

REF 01140	Product container filter with indicator	Size 190 mm
-----------	--	-----------------------

Product structure: sterilization crepe paper white, 60 g/m² according to DIN EN ISO 11607-1

Performance: round with diameter = 190 mm, with tab

Indicator: steam indicator Type 1 according to DIN EN ISO 11140-1:2014
colour change from pink to dark brown - black

Printing colour: 1 – colour, black, sterilizable

Toxicity: toxic effect is not known

Product description: Container filters made of paper for single use
for Aesculap-Containers
Round filters with steam indicator and a perforated hole
The Container filters are for single use and not recyclable.

Instruction for use: Place one filter over each perforated section on the inside of the SterilContainer lid. Secure each filter with the filter retention plate.

Compatibility

reusable sterilization containers:

Aesculap: ☒ suitable ☐ not suitable

Wagner: ☐ suitable ☒ not suitable

Martin: ☐ suitable ☒ not suitable


Medicon: ☐ suitable ☒ not suitable

SPECIFICATION

Sterilisation prozess:

STEAM	<input checked="" type="checkbox"/> suitable	<input type="checkbox"/> not suitable
EO	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable
FORM	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable
VH2O2	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> not suitable

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May 6th, 2021**DECLARATION OF CONFORMITY**

To whom it may concern,

We, Terragene S.A., located in Ruta Nacional N° 9, Km 280, Parque Industrial Micropi, Alvear (C.P. 2130), Santa Fe, Argentina, hereby declare that:

- Terragene Bionova® family of Biological Indicators for monitoring sterilization processes are developed, tested and produced in conformity with the following standards: ISO 11138-1:2017, ISO 11138-2:2017, ISO 11138-3:2017, ISO 11138-4:2017 and ISO 11138-5:2017. Chemical Indicator printed on SCBIs` labels are compliance with ISO 11140-1:2014 standard.
- Terragene Chemdye®, Integron® and Cintape® family of Chemical Indicators for monitoring sterilization processes are developed, tested and produced in conformity with the following standards: ISO 11140-1:2014, ISO 11140-3:2007, ISO 11140-4:2007 and ISO 11140-5:2007.
- Terragene Chemdye®, Integron® and Chemink® family of Cleaning and Hygiene Indicators are developed, tested and produced in conformity with the following standards: ISO 15883-1:2006 and ISO/TS 15883-5:2005.

The Biological and Chemical Indicators mentioned above are classified by ANMAT (National Administration for Drugs, Food and Medical Devices) disposition number 2318/2002 (to 2004) as Class 1 products according to the inherent risk posed to consumer's health. Cleaning and Hygiene Indicators are not considered medical devices in Argentina.

Our products are manufactured under ISO 13485:2016 / NS-EN ISO 13485:2016 Quality Management System Standards with certification number 243953-2017-AQ-ARG-NA-PS, granted by DNV international certification body.

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Magdalena Łochowska-Piotrowicz

splat CDWA4

**Trazanto****Cleaning indicator***For cleaning performance tests.*

Usage

Designed for routine control of cleaning processes and for monitoring automatic washing machines effectiveness as well as monitoring ultrasonic washing machines which operate at frequencies equal to or greater than 35 kHz, and temperatures between 30 °C and 70 °C.

Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.
ISO/TS 15883-5:2005.

Characteristics

Challenge level: very high.

Allows detection of all factors affecting the outcome of cleaning process.

Initial Color: fuchsia.

Size: 70 mm x 30 mm (approximate).

Synthetic temperature-stable support.

Result: elimination of fuchsia organic substance.

To interpret results visually, use the results reference guide on product's instructions for use.

Optionally, the readout of the cleaning indicator can be done using Trazanto Automatic Quality Control and Traceability System (see instructions for use in device's manual).

NOTE: If unfavorable results are obtained, check quality/quantity of detergent, cleaning time and temperature, positioning of the load, cleaning indicator location and general condition of the washing machine.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

Storage conditions

T = 10-30 °C, RH = 30-80 %, store in dark place within its original package.

Transportation conditions

Storage conditions should be strictly followed.

Transport in closed and strengthened boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

5 years.

Packing

200 units per aluminum foil bag.

Packing information: product description, intended use, instructions for use, standards reference, storage conditions, manufacturer information.

Labelling

On product's label: product's code and description, product's presentation, batch number, manufacture and expiration date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information

Read product's instructions for use thoroughly before use.

Description of identified and/or predictable risks which could represent a risk for the product

Store the product within its original package, away from light, disinfectants, gases and humidity.

Do not re-use the product.

Do not touch indicator ink with your hands.

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Thermodisinfection and cleaning indicator

For moist heat disinfection processes in washer-disinfectors

EN

Product description

After cleaning, surgical instruments must be disinfected to ensure safe subsequent handling and processing. Thermal disinfection with moist heat, is the most common method for disinfection of medical devices in the hospital setting. Disinfection effects can be defined by the A_0 value, a physical parameter specified in ISO 15883-1, which denotes the inactivation of microorganisms. A particular A_0 value can be achieved with the most diverse temperature/time combinations. For processing surgical instruments, an A_0 value of 600 is specified in ISO 15883-2, corresponding to a hold time of 1 minute at 90 °C, ensuring lethality of bacteria, mycobacteria, fungi and heat-sensitive viruses. If efficacy is also to be ensured against heat-resistant viruses, e.g. hepatitis B, many authorities in the field of disinfection, like the Robert Koch Institute, recommend an A_0 value of 3000, corresponding to a temperature of 90 °C with a hold time of 5 min. A higher safety approach from this institute comprises a cycle at 93°C, 10 minutes; a program known as the "BGA/RKI" or epidemic program. Integron® indicators have been designed to react in moist heat disinfection processes in washer-disinfectors, providing an accurate, convenient method of routine control, and assuring disinfection efficiency. The green indicating ink was developed to turn to purple when temperature and time conditions are met.

Characteristics

Integron® indicators consist of one-sided-laminated strips of 30 mm x 70 mm printed with green reactive ink, designed with state-of-the-art technology to withstand conditions present during moist heat disinfection processes

100% Toxic Heavy Metals free

Indications for use

Use the correspondent Integron® Chemical Indicator for monitoring thermal disinfection processes in washer-disinfectors according to the following programs:

IT27W-1: 90°C, 1 minute $A_0 = 600$
IT27W-5: 90°C, 5 minutes $A_0 = 3000$
IT27W-10: 93°C, 10 minutes. BGA/RKI/Epidemic Program

Precautions

Do not expose this product to Steam, EO, Dry Heat, or any sterilization or disinfection processes other than Moist Heat Disinfection Processes in washer-disinfectors. Do not reuse.

Instructions for use

- 1- Place the indicator on the holder.
 - 2- Place the system on the tray along with the material.
 - 3- Run the selected disinfection program and once finished remove the indicator.
 - 4- Visually analyze indicators and record results. The indicator must change to the reference color for considering that the desired disinfection conditions were met. If unsatisfactory results are obtained, perform the actions indicated therein and repeat the disinfection process.
- Use one Integron indicator for each cycle.

Interpretation of results

The color change of the indicating ink from green to purple indicates that the thermal disinfection conditions were attained, thus allowing a quick and easy assessment of the success of the process. If no color change is observed, the load cannot be considered disinfected. Failures during a disinfection cycle can be a consequence of inadequate low temperature or cycle duration, which lowers its efficiency. It is not recommended to release the load until a satisfactory result is obtained.

Storage and shelf life

Indicators should be stored in dry places, away from light, at temperatures between 0°C and 4 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.

Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date as well as the batch of the product are specified on the package label and on each unit of product. Do not use indicators after their expiration date.

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Integron©
Tıbbi Cihaz ve Sağlık Hizmetleri Ltd. Şti.