



3.0

## Disposable Biopsy Forceps (for enteroscope,O.D.3.0mm)

Model Number	O.D.of Jaw(mm)	Length(mm)	Minimum endoscope channel	Remarks
HM /101-A-30/2300	3.0	2300	3.3mm/3.6mm endoscope channel	for enteroscope, uncoated, oval cup, without spike
HM /101-A-30/2300S	3.0	2300	3.3mm/3.6mm endoscope channel	for enteroscope, coated, oval cup, without spike
HM/101-B-30/2300	3.0	2300	3.3mm/3.6mm endoscope channel	for enteroscope, uncoated, oval cup, with spike
HM/101-B-30/2300S	3.0	2300	3.3mm/3.6mm endoscope channel	for enteroscope, coated, oval cup, with spike

## Jednorazowe kleszcze biopsyjne (śr. 3.0 mm)

Model	Śr. Zewn	Długość (mm)	Min. kanał endoskopu	Uwagi
HM/101-A-30/2300	3,0	2300	Kanał endoskopu 3,3mm/3,6mm	Do enteroskopu, niepowlekane, owalne łyżeczki, bez igły
HM/101-A-30/2300S	3,0	2300	Kanał endoskopu 3,3mm/3,6mm	Do enteroskopu, powlekane, owalne łyżeczki, bez igły
HM/101-B-30/2300	3,0	2300	Kanał endoskopu 3,3mm/3,6mm	Do enteroskopu, niepowlekane, owalne łyżeczki, z igłą
HM/101-B-30/2300S	3,0	2300	Kanał endoskopu 3,3mm/3,6mm	Do enteroskopu, powlekane, owalne łyżeczki, z igłą