

ventimyst<sup>®</sup>





ventimyst<sup>®</sup>





**Zaawansowane aktywne nawilżanie VentiMyst®  
zapewnia maksymalną wilgotność i zarządza nadmiernym skraplaniem**

**Stosowanie nawilżaczy oddechowych i obwodów oddechowych z podgrzewanym drenem do aktywnego nawilżania zapewnia wiele korzyści, ale może również prowadzić do problemów spowodowanych skraplaniem.**


**VentiMyst® rozwiązuje te problemy, zapobiegając gromadzeniu się skroplin w obwodzie oddechowych.**

**Korzyści ze zmniejszonego do minimum skraplania dzięki VentiMyst®**

- Mniejsze ryzyko zanieczyszczenia docierającego do respiratora lub pacjenta**
- Minimalizuje potrzebę otwierania lub przerywania obwodu w celu usunięcia nadmiernego opadu deszczu, zmniejszając w ten sposób ryzyko zapalenia płuc nabytego przez respirator**
- Mniejsze ryzyko zwiększonego oporu na przepływ gazu w wyniku gromadzenia się płynu w odnodze wydechowej lub filtrze**

## VentiMyst® Moisture Management Technology


Reduced build-up of rain-out compared with standard corrugated heated wire limbs.



Standard Corrugated Heated Wire Limb  
with rain-out pooling in the limb corrugations

ventiMyst®  
No build-up of rain-out







## ventiMyst® Adult

038-31-252

038-31-252	<b>VentiMyst® Adult Dual Heated Wire with Pressure Line - Kit 1</b> Adult 22mm Dual Heated Wire Breathing System 1.6m, VentiMyst® Corrugated Expiratory Limb, Y-Piece with MDI Port and Pressure Port, 0.5m Limb, Pressure Line, 22M/22M,15F Connector, 22F/22M Elbow & Auto-Fill Chamber	10
038-31-253	<b>VentiMyst® Adult Dual Heated Wire without Pressure Line - Kit 2</b> As above without Pressure Line	10



038-33-252



038-33-252/S

## ventiMyst® Neonatal

038-33-252	<b>VentiMyst® Neonatal Dual Heated Wire with Auto-Fill Chamber - Kit 1</b> Neonatal 15mm Dual Heated Wire Breathing System 1.6m, VentiMyst® Corrugated Expiratory Limb, Y-Piece with Suction Port, 0.3m Incubator Extension 15mm, 8M/8F Ported T-Pieces, 22M/22M,15F Connector, 8M/15M Connector, Adapter Tube, 22F/22M Elbow, 9-11F/22M,15F Connectors, 10M,8F/15M,11F Connectors, 15F/22F Heated Wire Adapter, Detachable Pressure Line, 0.15m Limb 9-11F/9-11F, 0.6m Limb 22F/22F, 0.6m Limb 15M,11F/22F & Auto-Fill Chamber	10
038-33-253	<b>VentiMyst® Neonatal Dual Heated Wire without Chamber - Kit 2</b> As above without Humidification Chamber	10
038-33-252/S	<b>VentiMyst® Neonatal Dual Heated Wire for SLE Ventilator with Auto-Fill Chamber - Kit 1</b> Neonatal 15mm Dual Heated Wire Breathing System 1.6m, VentiMyst® Corrugated Expiratory Limb, Y-Piece with Suction Port & Fixed Pressure Line, 0.3m Incubator Extension 10mm, 22M/22M,15F Connector, 8M/8F Ported T-Piece, 22F/22M Elbow, 15M-22M/15F Connectors, Red & Amber Flow Restrictors, 0.15m Limb 9-11F/9-11F, 0.6m Limb 9-11F/22F, 0.6m Limb 15M,11F/22F & Auto-Fill Chamber	10
038-33-253/S	<b>VentiMyst® Neonatal Dual Heated Wire for SLE Ventilator without Chamber - Kit 2</b> As above without Humidification Chamber	10

SLE is a registered trademark of SLE Ltd.