



Certificate of Analysis

1.11374.0100 Buffer tablets pH 6.8 for preparing buffer solution acc. to WEISE for the staining of blood smears

Batch FN1557574

	Spec. Values		Batch Values	
pH-value (1 tablet in 1 l of water)	6.75 - 6.95		6.83	
Crushing strength	≥ 50	N	≥ 50	N
Disintegration time	≤ 20	min	5	min
Suitability for blood-smear staining	conforms		conforms	

Date of release (DD.MM.YYYY) 04.02.2025

Expiry date (DD.MM.YYYY) 31.10.2028

Dr. Timm Schlegelmilch

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.